



## Live-Image-Kit

TV-Set for Optical Measuring Instruments

# Live-Image-Kit

The classical application for the Live-Image-Kit is the display of the eyepiece image of Möller-Wedel Optical autocollimators, testing telescopes, alignment telescopes or autocollimators on a screen.

## Hardware und software included in the scope of delivery

- Video-Mini PC (Windows 11) with mouse and keyboard
- Monitor
- C-mount adapter and eyepiece hull
- USB CMOS camera (monochrome)
- Software LiveViewer Lite

## Video-Mini-PC

The live image kit includes a fully-fledged mini PC. This can also be used alternatively (without CMOS camera and C -mount adapter) as a space-saving solution for running software such as ELCOdirect and INTOMATIK-S in combination with a corresponding measuring system from Möller-Wedel Optical GmbH.

## Image field size with different eyepieces

When using the Live-Image-Kit with one of the three possible Möller-Wedel Optical eyepieces arises due to the focal length of the C-mounts and the used eyepiece a different magnification ratio. The resulting image field size is bordered by CMOS-Chip with 1:1,8" of the camera. The following images show the size of the reticle in relation to the image section.



Image section with eyepiece  $f=10\text{ mm}$



Image section with eyepiece  $f=14,7\text{ mm}$



Image section with eyepiece  $f=25\text{ mm}$

Art.-No.	Description	Comment
229 931	Live-Image-Kit	
Following components can be ordered seperately		
229 933	Video-Mini-PC	(Windows11)
133 401 07	C-mount adapter	Reccomended in combination with eyepiece hull
280 385	eyepiece hull	Reccomended in combination with C-Mount
773 000 33	USB CMOS camera	Monochrome camera with 1:1.8"
229 934	Software LiveViewer Lite	

Möller-Wedel Optical GmbH

Rosengarten 10  
D-22880 Wedel

Tel.: +49 - 41 03 - 9 37 76 10

Fax: +49 - 41 03 - 9 37 76 60

[www.moeller-wedel-optical.com](http://www.moeller-wedel-optical.com)

e-mail: [info@moeller-wedel-optical.com](mailto:info@moeller-wedel-optical.com)